

**Safe-T-Vue® • Product No. 7201**

**QUALITY ASSURANCE LOT PERFORMANCE REPORT**  
 Safe-T-Vue® X Non-Reversible Temperature Indicator Product No. 7201

Lot # **S10 040214**  
 Expiration Date: April-2016

**Non - Chemistry**

Mark  for "Accepted" and X for "Not Accepted."

| Box                                            | Adh. Label                                | Foil Lid                                       | Close Function                                    | Fill Level                                     |
|------------------------------------------------|-------------------------------------------|------------------------------------------------|---------------------------------------------------|------------------------------------------------|
| <input checked="" type="checkbox"/> Copy       | <input checked="" type="checkbox"/> Copy  | <input checked="" type="checkbox"/> Copy       | <input checked="" type="checkbox"/> Formed        | <input checked="" type="checkbox"/> Acceptable |
| <input checked="" type="checkbox"/> Tack       | <input checked="" type="checkbox"/> Tack  | <input checked="" type="checkbox"/> Position   | <input checked="" type="checkbox"/> Locked/Stable | <input type="checkbox"/> Underfill             |
| <input checked="" type="checkbox"/> Clean      | <input checked="" type="checkbox"/> Color | <input checked="" type="checkbox"/> Clean Peel | <b>Hinge</b>                                      | <input type="checkbox"/> Overfill              |
| <input checked="" type="checkbox"/> Box Damage | <b>Tray</b>                               | <b>Disk</b>                                    | <input checked="" type="checkbox"/> Elasticity    | <b>Color Label</b>                             |
|                                                | <input checked="" type="checkbox"/> Edges | <input checked="" type="checkbox"/> Appearance |                                                   | <input checked="" type="checkbox"/> Color      |
|                                                | <input checked="" type="checkbox"/> Color | <input checked="" type="checkbox"/> Fit        |                                                   | <input checked="" type="checkbox"/>            |

**Test Protocol – Finished Goods Quality Assurance**

Test performed on random sampling of finished product per ANSI ASQZ1.4-2003, General Inspection Level 1, Double Sampling Plan for Attributes with AQL=1%. Electronic thermometers with probes are calibrated annually. Download Validation SOP, from [www.williamlabs.com](http://www.williamlabs.com).

| <b>Temperature Change and Corresponding Safe-T-Vue® Change</b> |                                                              |
|----------------------------------------------------------------|--------------------------------------------------------------|
| Product Specification                                          | Safe-T-Vue 6.0±0.3°C, Safe-T-Vue 8±0.3°, Safe-T-Vue 10±0.3°C |
| Starting Test Temperatures                                     | 1 to 3 °C                                                    |

| Test Parameter                  | AQL | Sample Size | Accept Number | Test Result | Pass/Fail |
|---------------------------------|-----|-------------|---------------|-------------|-----------|
| Full Color Response Temperature | 1%  | 125         | 46150         | Pass        | Pass      |

**Performance Data:** Red Spots - Random 9.3 °C; Red Filling in White Center 9.7 °C; Full Red Color at 10.0 °C.

QUALITY ASSURANCE: APPROVED LOT # S10040214;

  
 Quality Manager

01May14  
 Date